



F198A mouse mAb

Catalog No	YP-mAb-09404
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	FAM198A C3orf41
Protein Name	Protein FAM198A
Immunogen	Synthesized peptide derived from human F198A AA range: 486-536
Specificity	This antibody detects endogenous levels of human F198A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted . Endoplasmic reticulum . Golgi apparatus . Membrane, caveola . Requires caveolae biogenesis to be secreted from the endoplasmic reticulum going through the Golgi apparatus where is post-translationally processed to the mature form. .
Tissue Specificity	Expressed in skin, lung and colon (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images