



CASK mouse mAb

Catalog No	YP-mAb-09406
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CSN3 CASK CSN10 CSNK
Protein Name	Kappa-casein
Immunogen	Synthesized peptide derived from human CASK AA range: 79-129
Specificity	This antibody detects endogenous levels of human CASK
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted.
Tissue Specificity	Mammary gland specific. Secreted in milk.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images