



RNF31 mouse mAb

Catalog No	YP-mAb-09419
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RNF31 ZIBRA
Protein Name	E3 ubiquitin-protein ligase RNF31 (EC 6.3.2.-) (HOIL-1-interacting protein) (HOIP) (RING finger protein 31) (Zinc in-between-RING-finger ubiquitin-associated domain protein)
Immunogen	Synthesized peptide derived from human RNF31 AA range: 833-883
Specificity	This antibody detects endogenous levels of human RNF31
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm .
Tissue Specificity	Expressed in both normal and transformed breast epithelial cell lines.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images