



CG023 mouse mAb

Catalog No	YP-mAb-09421
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C7orf23
Protein Name	Transmembrane protein C7orf23
Immunogen	Synthesized peptide derived from human CG023 AA range: 28-78
Specificity	This antibody detects endogenous levels of human CG023
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images