



ZBED1 mouse mAb

Catalog No	YP-mAb-09424
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZBED1 ALTE DREF KIAA0785 TRAMP
Protein Name	Zinc finger BED domain-containing protein 1 (Putative Ac-like transposable element) (dREF homolog)
Immunogen	Synthesized peptide derived from human ZBED1 AA range: 329-379
Specificity	This antibody detects endogenous levels of human ZBED1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus, PML body . Nucleus . In granular structures. .; Nucleus, PML body. (Microbial infection) Upon interaction with human adenovirus early E1A protein, the protein is redistributed to the peripheral areas of PML bodies. .
Tissue Specificity	Ubiquitously expressed at low levels (PubMed:9887332). Expression is highest in skeletal muscle, heart, spleen and placenta (PubMed:9887332).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images