



ARPC2 mouse mAb

Catalog No	YP-mAb-09425
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ARPC2 ARC34 PRO2446
Protein Name	Actin-related protein 2/3 complex subunit 2 (Arp2/3 complex 34 kDa subunit) (p34-ARC)
Immunogen	Synthesized peptide derived from human ARPC2 AA range: 147-197
Specificity	This antibody detects endogenous levels of human ARPC2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm, cytoskeleton . Cell projection . Cell junction, synapse, synaptosome . Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images