



GEMI4 mouse mAb

Catalog No	YP-mAb-09435
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GEMIN4
Protein Name	Gem-associated protein 4 (Gemin-4) (Component of gems 4) (p97)
Immunogen	Synthesized peptide derived from human GEMI4 AA range: 754-804
Specificity	This antibody detects endogenous levels of human GEMI4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm. Nucleus. Nucleus, nucleolus. Nucleus, gem. Localized in subnuclear structures next to coiled bodies, called gems, which are highly enriched in spliceosomal snRNPs and in the nucleolus.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images