



PALM3 mouse mAb

Catalog No	YP-mAb-09489
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PALM3
Protein Name	Paralemmin-3
Immunogen	Synthesized peptide derived from human PALM3 AA range: 132-182
Specificity	This antibody detects endogenous levels of human PALM3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Cell membrane ; Lipid-anchor .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images