



REEP6 mouse mAb

Catalog No	YP-mAb-09502
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	REEP6 C19orf32 DP1L1
Protein Name	Receptor expression-enhancing protein 6 (Polyposis locus protein 1-like 1)
Immunogen	Synthesized peptide derived from human REEP6 AA range: 110-160
Specificity	This antibody detects endogenous levels of human REEP6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Cytoplasmic vesicle, clathrin-coated vesicle membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in circumvallate papillae and testis (PubMed:16720576). Expressed in the retina. Isoform 1 is predominantly present in mature optic cups. Isoform 1 expression is confined to the cell body and inner segment of developing rod photoreceptor cells (PubMed:27889058).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images