



NXF5 mouse mAb

Catalog No	YP-mAb-09508
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	NXF5 TAPL1
Protein Name	Nuclear RNA export factor 5 (TAP-like protein 1) (TAPL-1)
Immunogen	Synthesized peptide derived from human NXF5 AA range: 245-295
Specificity	This antibody detects endogenous levels of human NXF5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm. Nucleus. Mainly localized in the cytoplasm of cells and more particularly in the cell body and neurites of hippocampal neurons. Although nuclear localization is also observed. Not detected at nuclear rim.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images