



NDK5 mouse mAb

Catalog No	YP-mAb-09531
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	NME5
Protein Name	Nucleoside diphosphate kinase homolog 5 (NDK-H 5) (NDP kinase homolog 5) (Inhibitor of p53-induced apoptosis-beta) (IPIA-beta) (Testis-specific nm23 homolog) (nm23-H5)
Immunogen	Synthesized peptide derived from human NDK5 AA range: 25-75
Specificity	This antibody detects endogenous levels of human NDK5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell projection, cilium .
Tissue Specificity	Specifically expressed in testis germinal cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images