



RDH13 mouse mAb

Catalog No	YP-mAb-09542
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RDH13 PSEC0082 UNQ736/PRO1430
Protein Name	Retinol dehydrogenase 13 (EC 1.1.1.-)
Immunogen	Synthesized peptide derived from human RDH13 AA range: 1-51
Specificity	This antibody detects endogenous levels of human RDH13
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Mitochondrion inner membrane ; Peripheral membrane protein . Localized on the outer side of the inner mitochondrial membrane. .
Tissue Specificity	Widely expressed (PubMed:18039331). In the retina, detected in the inner segment of the photoreceptor cells. Weak signals are observed in a small population of inner nuclear neurons and the inner plexiform layer (PubMed:12226107).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images