



BTBD8 mouse mAb

Catalog No	YP-mAb-09546
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	BTBD8
Protein Name	BTB/POZ domain-containing protein 8
Immunogen	Synthesized peptide derived from human BTBD8 AA range: 224-274
Specificity	This antibody detects endogenous levels of human BTBD8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Localized to nucleus in fetal cells.
Tissue Specificity	Highly expressed in fetal brain. Weakly expressed in adult brain and prostate.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images