



TIM23 mouse mAb

Catalog No	YP-mAb-09562
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TIMM23 TIM23
Protein Name	Mitochondrial import inner membrane translocase subunit Tim23
Immunogen	Synthesized peptide derived from human TIM23 AA range: 95-145
Specificity	This antibody detects endogenous levels of human TIM23
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	23kD
Cell Pathway	Mitochondrion inner membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images