



NUMBL mouse mAb

Catalog No	YP-mAb-09564
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NUMBL
Protein Name	Numb-like protein (Numb-related protein) (Numb-R)
Immunogen	Synthesized peptide derived from human NUMBL AA range: 290-340
Specificity	This antibody detects endogenous levels of human NUMBL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Symmetrically distributed throughout the cytoplasm in non dividing neuroblasts of the CNS. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images