



VAMP5 mouse mAb

Catalog No	YP-mAb-09580
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VAMP5 HSPC191
Protein Name	Vesicle-associated membrane protein 5 (VAMP-5) (Myobrevin)
Immunogen	Synthesized peptide derived from human VAMP5 AA range: 2-52
Specificity	This antibody detects endogenous levels of human VAMP5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell membrane ; Single-pass type IV membrane protein . Endomembrane system ; Single-pass type IV membrane protein . Golgi apparatus, trans-Golgi network membrane ; Single-pass type IV membrane protein . Associated with the plasma membrane as well as intracellular perinuclear and peripheral vesicular structures of myotubes. Associated with the trans-Golgi, but not with the cis-Golgi apparatus (By similarity). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images