



CT191 mouse mAb

Catalog No	YP-mAb-09583
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	NCOR1P1 C20orf191
Protein Name	Putative nuclear receptor corepressor 1-like protein NCOR1P1 (Nuclear receptor corepressor 1 pseudogene 1)
Immunogen	Synthesized peptide derived from human CT191 AA range: 12-62
Specificity	This antibody detects endogenous levels of human CT191
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images