



ZN644 mouse mAb

Catalog No	YP-mAb-09600
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF644 KIAA1221 ZEP2
Protein Name	Zinc finger protein 644 (Zinc finger motif enhancer-binding protein 2) (Zep-2)
Immunogen	Synthesized peptide derived from human ZN644 AA range: 674-724
Specificity	This antibody detects endogenous levels of human ZN644
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	146kD
Cell Pathway	Nucleus.
Tissue Specificity	Expressed in liver, placenta, retina and retinal pigment epithelium.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images