



SCC4 mouse mAb

Catalog No	YP-mAb-09609
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MAU2 KIAA0892 SCC4
Protein Name	MAU2 chromatid cohesion factor homolog (MAU-2) (Cohesin loading complex subunit SCC4 homolog)
Immunogen	Synthesized peptide derived from human SCC4 AA range: 558-608
Specificity	This antibody detects endogenous levels of human SCC4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus, nucleoplasm . Nucleus . Chromosome . Binds to chromatin from the end of mitosis until prophase.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images