



REEP2 mouse mAb

Catalog No	YP-mAb-09611
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	REEP2 C5orf19 SGC32445
Protein Name	Receptor expression-enhancing protein 2
Immunogen	Synthesized peptide derived from human REEP2 AA range: 195-245
Specificity	This antibody detects endogenous levels of human REEP2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Detected in brain, heart and skeletal muscle, and at low levels in placenta, kidney and pancreas (PubMed:11161817). Expressed in circumvallate papillae (PubMed:16720576).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images