



K1210 mouse mAb

Catalog No	YP-mAb-09642
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	KIAA1210
Protein Name	Uncharacterized protein KIAA1210
Immunogen	Synthesized peptide derived from human K1210 AA range: 747-797
Specificity	This antibody detects endogenous levels of human K1210
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasmic vesicle, secretory vesicle, acrosome . Localizes to the sex body in spermatocyte, acrosome, and near the ectoplasmic specialization. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images