



TPGS1 mouse mAb

Catalog No	YP-mAb-09644
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TPGS1 C19orf20
Protein Name	Tubulin polyglutamylase complex subunit 1 (PGs1)
Immunogen	Synthesized peptide derived from human TPGS1 AA range: 158-208
Specificity	This antibody detects endogenous levels of human TPGS1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm, cytoskeleton, cilium axoneme . Cytoplasm, cytoskeleton, flagellum axoneme . Cytoplasm, cytoskeleton, cilium basal body . Cytoplasm, cytoskeleton, flagellum basal body . Cell projection, axon . Cell projection, dendrite . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriolar satellite . Associated with microtubules from neurites, centrosomes, basal bodies and axonemes. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images