



D107A mouse mAb

Catalog No	YP-mAb-09670
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	DEFB107A DEFB107 DEFB7; DEFB107B
Protein Name	Beta-defensin 107 (Beta-defensin 7) (BD-7) (DEFB-7) (Defensin, beta 107)
Immunogen	Synthesized peptide derived from human D107A AA range: 3-53
Specificity	This antibody detects endogenous levels of human D107A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted.
Tissue Specificity	Specifically expressed in testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images