



SAC31 mouse mAb

Catalog No	YP-mAb-09676
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SAC3D1 SHD1
Protein Name	SAC3 domain-containing protein 1 (SAC3 homology domain-containing protein 1)
Immunogen	Synthesized peptide derived from human SAC31 AA range: 31-81
Specificity	This antibody detects endogenous levels of human SAC31
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, spindle . Localizes on centrosomes in interphase cells and at spindles in mitosis. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images