



POK7 mouse mAb

Catalog No	YP-mAb-09685
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	
Protein Name	HERV-K_1q23.3 provirus ancestral Pol protein (HERV-K(C1a) Pol protein) (HERV-K110 Pol protein) (HERV-K18 Pol protein) [Includes: Reverse transcriptase (EC 2.7.7.49); Ribonuclease H (RNase H) (EC 3.1.2
Immunogen	Synthesized peptide derived from human POK7 AA range: 450-500
Specificity	This antibody detects endogenous levels of human POK7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images