



CPT2 mouse mAb

Catalog No	YP-mAb-09692
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CPT2 CPT1
Protein Name	Carnitine O-palmitoyltransferase 2, mitochondrial (EC 2.3.1.21) (Carnitine palmitoyltransferase II) (CPT II)
Immunogen	Synthesized peptide derived from human CPT2 AA range: 459-509
Specificity	This antibody detects endogenous levels of human CPT2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	72kD
Cell Pathway	Mitochondrion inner membrane; Peripheral membrane protein; Matrix side.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images