



PRAF1 mouse mAb

Catalog No	YP-mAb-09696
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RABAC1 PRA1 PRAF1
Protein Name	Prenylated Rab acceptor protein 1 (PRA1 family protein 1)
Immunogen	Synthesized peptide derived from human PRAF1 AA range: 130-180
Specificity	This antibody detects endogenous levels of human PRAF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Cytoplasm . Golgi apparatus . Cytoplasmic vesicle, secretory vesicle, synaptic vesicle . According to some authors, it is an integral membrane protein, while others showed that it is cytoplasmic and membrane-associated to Golgi and synaptic vesicles. .
Tissue Specificity	Ubiquitous. Strongest expression found in placenta, pituitary gland, kidney, lung and stomach.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images