



CYTSB mouse mAb

Catalog No	YP-mAb-09709
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SPECC1 CYTSB NSP5
Protein Name	Cytospin-B (Nuclear structure protein 5) (NSP5) (Sperm antigen HCMOGT-1) (Sperm antigen with calponin homology and coiled-coil domains 1)
Immunogen	Synthesized peptide derived from human CYTSB AA range: 403-453
Specificity	This antibody detects endogenous levels of human CYTSB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .; [Isoform 3]: Membrane ; Lipid-anchor .; [Isoform 4]: Membrane ; Lipid-anchor .; [Isoform 5]: Membrane ; Lipid-anchor .
Tissue Specificity	Highly expressed in testis. Barely detectable in other tissues. Also highly expressed in some cancer cell lines.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images