



ODFP4 mouse mAb

Catalog No	YP-mAb-09710
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ODF4 OPPO1
Protein Name	Outer dense fiber protein 4 (Outer dense fiber of sperm tails protein 4) (Testis-specific protein oppo 1) (hOPPO1)
Immunogen	Synthesized peptide derived from human ODFP4 AA range: 119-169
Specificity	This antibody detects endogenous levels of human ODFP4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in testis and sperm; especially localized to sperm tail (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images