



GALP mouse mAb

Catalog No	YP-mAb-09719
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GALP
Protein Name	Galanin-like peptide
Immunogen	Synthesized peptide derived from human GALP AA range: 12-62
Specificity	This antibody detects endogenous levels of human GALP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted.
Tissue Specificity	Isoform 2 is found in ganglia of ganglioneuroma and ganglioneuroblastoma, as well as in differentiated tumor cells of neuroblastoma tissues. Not found in undifferentiated neuroblasts. Isoform 2 is found in the skin, in pericytes covering microvascular arterioles and venules on their abluminal surfaces. In larger vessels, isoform 2 is expressed in layers of smooth muscle cells. Isoform 2 is not detected in endothelial cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images