



TRH mouse mAb

Catalog No	YP-mAb-09721
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRH
Protein Name	Prothyroliberin [Cleaved into: Thyroliberin (Protirelin) (TSH-releasing factor) (Thyrotropin-releasing factor) (TRF) (Thyrotropin-releasing hormone) (TRH)]
Immunogen	Synthesized peptide derived from human TRH AA range: 133-183
Specificity	This antibody detects endogenous levels of human TRH
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted.
Tissue Specificity	Hypothalamus. Expressed in the hair follicle epithelium (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images