



UD2B7 mouse mAb

Catalog No	YP-mAb-09725
Isotype	IgG
Reactivity	Human;Rat
Applications	WB
Gene Name	UGT2B7 UGTB2B9
Protein Name	UDP-glucuronosyltransferase 2B7 (UDPGT 2B7) (EC 2.4.1.17) (3,4-catechol estrogen-specific UDPGT) (UDP-glucuronosyltransferase 2B9) (UDPGT 2B9) (UDPGTh-2)
Immunogen	Synthesized peptide derived from human UD2B7 AA range: 467-517
Specificity	This antibody detects endogenous levels of human UD2B7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images