



BDP1 mouse mAb

Catalog No	YP-mAb-09727
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	BDP1 KIAA1241 KIAA1689 TFNR
Protein Name	Transcription factor TFIIIB component B" homolog (Transcription factor IIIB 150) (TFIIIB150) (Transcription factor-like nuclear regulator)
Immunogen	Synthesized peptide derived from human BDP1 AA range: 1732-1782
Specificity	This antibody detects endogenous levels of human BDP1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	Isoform 2 is highly expressed in cerebellum.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images