



COX7S mouse mAb

Catalog No	YP-mAb-09732
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	COX7A2P2 COX7A3 COX7AL2 COX7AP2
Protein Name	Putative cytochrome c oxidase subunit 7A3, mitochondrial (Cytochrome c oxidase subunit VIIa 3)
Immunogen	Synthesized peptide derived from human COX7S AA range: 13-63
Specificity	This antibody detects endogenous levels of human COX7S
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Mitochondrion inner membrane .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images