



CLU mouse mAb

Catalog No	YP-mAb-09738
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	KIAA0664
Protein Name	Clustered mitochondria protein homolog
Immunogen	Synthesized peptide derived from human CLU AA range: 435-485
Specificity	This antibody detects endogenous levels of human CLU
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Cytoplasmic granule . A fraction colocalizes with tyrosinated tubulin and can be detected close to mitochondria. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images