



UBR4 mouse mAb

Catalog No	YP-mAb-09764
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	UBR4 KIAA0462 KIAA1307 RBAF600 ZUBR1
Protein Name	E3 ubiquitin-protein ligase UBR4 (EC 6.3.2.-) (600 kDa retinoblastoma protein-associated factor) (N-recognin-4) (Retinoblastoma-associated factor of 600 kDa) (RBAF600) (p600) (Zinc finger UBR1-type pr
Immunogen	Synthesized peptide derived from human UBR4 AA range: 809-859
Specificity	This antibody detects endogenous levels of human UBR4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Multi-pass membrane protein . Cytoplasm. Cytoplasm, cytoskeleton. Nucleus. Concentrates at the leading edge of membrane structures involved in actin motility.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images