



MCA3 mouse mAb

Catalog No	YP-mAb-09770
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	EEF1E1 AIMP3 P18
Protein Name	Eukaryotic translation elongation factor 1 epsilon-1 (Aminoacyl tRNA synthetase complex-interacting multifunctional protein 3) (Elongation factor p18) (Multisynthase complex auxiliary component p18)
Immunogen	Synthesized peptide derived from human MCA3 AA range: 70-120
Specificity	This antibody detects endogenous levels of human MCA3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Cytoplasm, cytosol . Nucleus . Cytoplasmic under growth arrest conditions. Translocated into the nucleus when growth resumes (S phase) and following DNA damage.
Tissue Specificity	Down-regulated in various cancer tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images