



MTA3 mouse mAb

Catalog No	YP-mAb-09790
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MTA3 KIAA1266
Protein Name	Metastasis-associated protein MTA3
Immunogen	Synthesized peptide derived from human MTA3 AA range: 270-320
Specificity	This antibody detects endogenous levels of human MTA3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in germinal centers of lymphoid tissues. No expression in nonepithelial cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images