



K2C80 mouse mAb

Catalog No	YP-mAb-09793
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KRT80 KB20
Protein Name	Keratin, type II cytoskeletal 80 (Cytokeratin-80) (CK-80) (Keratin-80) (K80) (Type-II keratin Kb20)
Immunogen	Synthesized peptide derived from human K2C80 AA range: 107-157
Specificity	This antibody detects endogenous levels of human K2C80
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	Weakly expressed in tongue, but not skin or in any other tissues or organs examined.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images