



ERI3 mouse mAb

Catalog No	YP-mAb-09811
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ERI3 PINT1 PRNPIP PRNPIP1
Protein Name	ERI1 exoribonuclease 3 (EC 3.1.-.-) (Prion interactor 1) (Prion protein-interacting protein)
Immunogen	Synthesized peptide derived from human ERI3 AA range: 64-114
Specificity	This antibody detects endogenous levels of human ERI3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images