



LMAN2 mouse mAb

Catalog No	YP-mAb-09884
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	LMAN2 C5orf8
Protein Name	Vesicular integral-membrane protein VIP36 (Glycoprotein GP36b) (Lectin mannose-binding 2) (Vesicular integral-membrane protein 36) (VIP36)
Immunogen	Synthesized peptide derived from human LMAN2 AA range: 227-277
Specificity	This antibody detects endogenous levels of human LMAN2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Endoplasmic reticulum-Golgi intermediate compartment membrane ; Single-pass type I membrane protein . Golgi apparatus membrane ; Single-pass membrane protein . Endoplasmic reticulum membrane ; Single-pass type I membrane protein .
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images