



NAT8 mouse mAb

Catalog No	YP-mAb-10001
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NAT8 CML1 GLA TSC501
Protein Name	Probable N-acetyltransferase 8 (EC 2.3.1.-) (Camello-like protein 1)
Immunogen	Synthesized peptide derived from human NAT8 AA range: 38-88
Specificity	This antibody detects endogenous levels of human NAT8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Endoplasmic reticulum-Golgi intermediate compartment membrane; Single-pass type II membrane protein. Endoplasmic reticulum membrane; Single-pass type II membrane protein.
Tissue Specificity	Preferentially expressed in liver and kidney. Also detected in brain (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images