



LBR mouse mAb

Catalog No	YP-mAb-10021
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LBR
Protein Name	Lamin-B receptor (Integral nuclear envelope inner membrane protein) (LMN2R)
Immunogen	Synthesized peptide derived from human LBR AA range: 411-461
Specificity	This antibody detects endogenous levels of human LBR
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus inner membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane . Cytoplasm . Nucleus . Nucleus; nuclear rim. .
Tissue Specificity	Expressed in the bone marrow, liver, heart, adrenal gland, lung, placenta and uterus (PubMed:16784888). Expressed in osteoclasts and osteoblast-like cells (PubMed:21327084).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images