



TIM16 mouse mAb

Catalog No	YP-mAb-10027
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PAM16 MAGMAS TIM16 TIMM16 CGI-136
Protein Name	Mitochondrial import inner membrane translocase subunit TIM16 (Mitochondria-associated granulocyte macrophage CSF-signaling molecule) (Presequence translocated-associated motor subunit PAM16)
Immunogen	Synthesized peptide derived from human TIM16 AA range: 46-96
Specificity	This antibody detects endogenous levels of human TIM16
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Mitochondrion inner membrane ; Peripheral membrane protein ; Matrix side .
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images