



CC154 mouse mAb

Catalog No	YP-mAb-10030
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CCDC154
Protein Name	Coiled-coil domain-containing protein 154
Immunogen	Synthesized peptide derived from human CC154 AA range: 154-204
Specificity	This antibody detects endogenous levels of human CC154
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Early endosome .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images