



PGTB2 mouse mAb

Catalog No	YP-mAb-10031
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RABGGTB GGTB
Protein Name	Geranylgeranyl transferase type-2 subunit beta (EC 2.5.1.60) (Geranylgeranyl transferase type II subunit beta) (GGTase-II-beta) (Rab geranyl-geranyltransferase subunit beta) (Rab GG transferase beta)
Immunogen	Synthesized peptide derived from human PGTB2 AA range: 64-114
Specificity	This antibody detects endogenous levels of human PGTB2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images