



## MSL2 mouse mAb

|                           |                                                                                                                                                   |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No                | YP-mAb-10036                                                                                                                                      |
| Isotype                   | IgG                                                                                                                                               |
| Reactivity                | Human;Mouse                                                                                                                                       |
| Applications              | WB                                                                                                                                                |
| Gene Name                 | MSL2 KIAA1585 MSL2L1 RNF184                                                                                                                       |
| Protein Name              | Male-specific lethal 2 homolog (MSL-2) (Male-specific lethal 2-like 1) (MSL2-like 1) (Male-specific lethal-2 homolog 1) (RING finger protein 184) |
| Immunogen                 | Synthesized peptide derived from human MSL2 AA range: 9-59                                                                                        |
| Specificity               | This antibody detects endogenous levels of human MSL2                                                                                             |
| Formulation               | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                           |
| Source                    | Monoclonal, Mouse,IgG                                                                                                                             |
| Purification              | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                              |
| Dilution                  | WB 1:500-1:2000                                                                                                                                   |
| Concentration             | 1 mg/ml                                                                                                                                           |
| Purity                    | ≥90%                                                                                                                                              |
| Storage Stability         | -20°C/1 year                                                                                                                                      |
| Synonyms                  |                                                                                                                                                   |
| Observed Band             |                                                                                                                                                   |
| Cell Pathway              |                                                                                                                                                   |
| Tissue Specificity        |                                                                                                                                                   |
| Function                  |                                                                                                                                                   |
| Background                |                                                                                                                                                   |
| matters needing attention | Avoid repeated freezing and thawing!                                                                                                              |
| Usage suggestions         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.                 |



## Products Images