



PGAP3 mouse mAb

Catalog No	YP-mAb-10045
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PGAP3 CAB2 PERLD1 UNQ546/PRO1100
Protein Name	Post-GPI attachment to proteins factor 3 (COS16 homolog) (hCOS16) (Gene coamplified with ERBB2 protein) (PER1-like domain-containing protein 1)
Immunogen	Synthesized peptide derived from human PGAP3 AA range: 252-302
Specificity	This antibody detects endogenous levels of human PGAP3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Golgi apparatus membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane ; Multi-pass membrane protein . Mainly localizes to Golgi apparatus. .
Tissue Specificity	Ubiquitously expressed, with highest levels in thyroid and placenta.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images