



TRI17 mouse mAb

Catalog No	YP-mAb-10055
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRIM17 RBCC RNF16 TERF
Protein Name	E3 ubiquitin-protein ligase TRIM17 (EC 6.3.2.-) (RING finger protein 16) (Testis RING finger protein) (Tripartite motif-containing protein 17)
Immunogen	Synthesized peptide derived from human TRI17 AA range: 95-145
Specificity	This antibody detects endogenous levels of human TRI17
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Lysosome .
Tissue Specificity	Almost exclusively in the testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images