



PGPIL mouse mAb

Catalog No	YP-mAb-10059
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PGPEP1L
Protein Name	Pyroglutamyl-peptidase 1-like protein (EC 3.4.19.-)
Immunogen	Synthesized peptide derived from human PGPIL AA range: 115-165
Specificity	This antibody detects endogenous levels of human PGPIL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images